



## ZETA PLUS® VR SERIES Viral Reduction Depth Filters

*VR Series For Enhanced Log Reduction of Viruses From Biological Fluids*



### Application

The removal and/or inactivation of contaminating viruses from biotherapeutics is a requisite for ensuring final product safety. Zeta Plus® VR cartridge depth filters remove significant levels of viruses from biological fluids. They provide validatable viral titer reduction, high flow rates, scalability, economy, disposability and ease-of-use in the biological manufacturing environment. The VR Series includes specific filter media recommendations for virus removal from blood plasma proteins and bioprocess derived cell culture fluids.

Virus Removal Applications for Zeta Plus® VR Media :

- Processing of human plasma-derived proteins
- Recombinant proteins from cell culture and fermentation
- Monoclonal antibodies
- Serum and other raw materials derived from animal sources

Zeta Plus® depth filtration offers an attractive and economical alternative to other available technologies for reduction of viral burden. Zeta Plus® VR Series depth filters can be used as a bolt-on processing step to existing in-process viral clearance methodologies, providing an added measure of end-product safety.

### Features

- Combined adsorption & physical entrapment virus removal mechanisms
- Quality Control tested for media adsorption capacity
- Self-contained, single-use disposable capsules & cartridge modules
- Range of cartridge sizes, sanitary design cartridge housings & disposable capsules
- Optimised cartridge designs for biologicals filtration
- Complete technical, regulatory & validation support
- Economical “bolt-on” separation technology as adjunct or prefilter to other viral clearance steps

### Benefits

- Validatable viral log reduction
- Provides added level of viral safety
- Assurance of filter media consistency & performance
- Single-use, disposable design eliminates validation concerns associated with reusable clearance options (chromatography column sanitization & regeneration)
- Eliminates concern over virus carryover from one purification cycle to the next
- Facilitates scaleup from small volume feasibility trials & validation studies to production scale operations
- Enable high flow processes
- Robust in-situ steam-sterilisable cartridges
- Drug Master Files simplify FDA submissions
- Regulatory Support Files detail filter media regulatory compliance & safety data
- Provides added level of validatable viral clearance
- Can replace or augment chromatography columns or function as prefilter to final virus removal membrane filter

Regulators (FDA) recommend validation clearance via a combination of removal (filtration, chromatography) and inactivation (temperature, pH, solvent) steps where multiple methods of virus clearance are required to assure product safety.

## Zeta Plus® VR Filter Selection

A range of cartridge sizes is available for varying batch size, so scale-up is easy and predictable from lab-scale, to pilot plant, to full production scale. Three different cartridge sizes (8, 12 and 16" diameter) make scale-up straightforward and predictable. Zeta Plus® cartridges are used in a flow-through mode, enabling the process stream to be passed through the filter cartridge while viruses are retained within the depth filter matrix. Small, low area disposable BioCap® capsules are available for laboratory scale and process development work and are ideal for scaled down viral validation studies.

## Zeta Plus® Filter Mechanisms of Virus Removal

Zeta Plus® filter media is a family of cellulosic depth filtration media designed to retain contaminants by both physical entrapment and adsorption. It is composed of high-area process filter aids embedded in a cellulose fibre depth matrix. During the manufacturing process, a cationic charge modifier is chemically bound to the matrix components, forming a permanent, interconnected, rigid depth filter with positively-charged electrokinetic capture sites. The resulting porous depth filter structure is a tortuous network of adsorptive flow channels capable of retaining contaminating viruses through a combination of adsorption and mechanical entrapment. Cumulative viral titer reduction may also be attainable through a staged two-step Zeta Plus® system.

Zeta Plus® depth filtration is a key component in biopharmaceutical production schemes for clarifying cell culture fluids (cell separation) and is utilized worldwide in the fractionation of human and animal plasma proteins.

## Selection Guide - Zeta Plus® VR Filter Media

A range of Zeta Plus® VR depth filter media is available for viral reduction from plasma proteins and bioprocess-derived (cell culture) proteins. Specific recommendations for each application are provided in Table 1 below. VR filter media is available in a range of porosities. Each VR filter media is tested and optimised for biologicals and bioprocess filtration (pharmaceutical-grade).

All cartridge and capsule components have been USP Class VI tested for biological safety.

**Table 1**

Selection of Zeta Plus® VR Series Filter Media

Filter Media Choice	Recommended for:		Nominal Retention Rating (µm)	20°C Water Flow @ 345 mb (Lpm/m²)	Auto-clavable or In-situ Steam Sterilisable	[API <sup>31</sup> ] in 20°C water @ 54 l/m²
	Plasma Proteins	Cell Culture Proteins				
VR 01	◆		0.8	50	No	40
VR 02	◆	◆	0.8	77	Yes	20
VR 03		◆	0.8	50	Yes	40
VR 04	◆		0.45	37	No	20
VR 05		◆	0.2	16	Yes	40
VR 06	◆	◆	0.2	24	Yes	20
VR 07	◆	◆	0.2	16	Yes	60

## Viral Reduction Validation

Many factors influence the choice of process steps to be studied when performing viral clearance evaluation and validation. Viral clearance validation studies are performed with a select panel of model viruses. Different virus panels are chosen to validate blood and plasma products, and separately, cell-culture derived products.

Virus removal efficiency may be dependent on fluid and processing conditions and is best determined during qualification and validation studies using suitable virus model systems in a controlled laboratory study. Contact CUNO Scientific Applications Support Services (SASS) for further technical advice concerning validation of viral reduction with Zeta Plus® VR depth filters.

Data is available which confirms the ability of Zeta Plus® VR depth filter media to function as an effective and validatable viral log reduction filter. Table 2 illustrates results for a reported viral clearance study utilising Zeta Plus® VR depth filters.

**Table 2**

Viral Clearance In An Immunoglobulin Production Process

Process Step	Cumulative Virus Titer Reduction (Log <sub>10</sub> )				
	BVD	EMC	HIV	PPV	PRV
Solvent detergent	> 4.3	—	> 5.3	—	> 7.3
Supernatant III	1.4	4.3	6.1	4.7	3.8
Zeta Plus® VR 03 Depth Filtration	4.8	4.5	4.7	3.7	5.4
Total Cumulative Reduction	> 10.5	8.8	> 16.1	8.4	> 16.6

Source : D. Revie, *Novel Validation Approaches to Obtain Maximum Viral Clearance from an Immunoglobulin Production Process*, IBC 2nd International Symposium on Viral Clearance, Philadelphia, PA, June, 1998.

## Quality Control of Zeta Plus® VR Series Filters

CUNO applies rigorous Quality Control testing and standards both during in-process manufacturing and during final lot release. This testing ensures consistent filter media performance in critical virus removal applications. A Certificate of Quality is provided with each VR Series filter cartridge and disposable BioCap® capsule.

Each grade of Zeta Plus® VR Series filter media is Quality Control tested on a lot release basis for the presence and magnitude of positive charge.

Summarized in Table 3 are key quality control test parameters and lot release test criteria applied to VR Series filters.

**Table 3**  
Quality Control of Zeta Plus® VR Series Filters

Filter media (µm)	* Test for Filter Media Charge Capacity	* Optimized & Tested for Low AI Extractables	* Filter Materials Optimized for Low LAL-RM	Non-pyrogenic per USP Endotoxin Test	100% 21 CFR Listed Materials	Supplied with Certificate of Quality	USP Class VI Biological Tests
VR 01 (0.8)	Yes	No	No	Yes	Yes	Yes	Pass
VR 02 (0.8)	Yes	Yes	Yes	Yes	Yes	Yes	Pass
VR 03 (0.8)	Yes	No	No	Yes	Yes	Yes	Pass
VR 04 (0.45)	Yes	Yes	No	Yes	Yes	Yes	Pass
VR 05 (0.2)	Yes	No	No	Yes	Yes	Yes	Pass
VR 06 (0.2)	Yes	Yes	Yes	Yes	Yes	Yes	Pass
VR 07 (0.2)	Yes	Yes	No	Yes	Yes	Yes	Pass

\* Indicates filter media manufacturing lot is sampled and tested as lot release test criteria

## Zeta Plus® VR Disposable Capsules, Filter Cartridges, and Filter Sheets - Configured to Suit Your Every Need

Zeta Plus® VR Series filter products are highly scaleable and can be used for applications ranging from millilitre lab volumes to production scale operations. Filter media is available in the following configurations:

### VR Filter Sheets

VR filter media is available in a range of filter sheet sizes to fit plate-and-frame filter presses (contact CUNO for specific sheet sizes and part numbers).

### BioCap® Disposable VR Filter Capsules

VR filter media is available in three different sizes of BioCap® disposable capsules for bench scale and pilot scale filtration (filter areas of 25 cm<sup>2</sup>, 650 cm<sup>2</sup>, and 1300 cm<sup>2</sup>). For more information on BioCap® disposable assemblies, including operating specifications.

### VR Filter Cartridges and Sanitary Housings

VR filter media is available in 8", 12" and 16" diameter, easy-to-use, disposable cartridges, providing surface areas from 650 cm<sup>2</sup> to 3.2 m<sup>2</sup> per cartridge. A range of sanitary design Zeta Plus® filter housings is available for VR Series cartridges.



Zeta Plus® VR Series filter media is available in a wide range of cartridge sizes including 8", 12", and 16" diameter filter cartridges for varying flow rates and batch volumes. CUNO's sanitary design ZPC and ZPB Zeta Plus® filter housings are designed for pharmaceutical and biological manufacturing environments. For more information on Zeta Plus® sanitary filter housings, contact CUNO.

## SCIENTIFIC APPLICATIONS SUPPORT SERVICES (SASS)

CUNO's 90+ years of experience are synonymous with quality, performance and high level technical support. The cornerstone of CUNO's philosophy is service to customers, not only in product quality and prompt delivery, but also in validation, in applications support and in the sharing of scientific information. CUNO's Scientific Applications Support Services group works closely with customers to solve difficult separation problems and to recommend the most economical and efficient filter systems. SASS specialists routinely work with third-party contract laboratories to validate Zeta Plus® filters for viral clearance applications. Contact CUNO's SASS group for more details on viral validation support.


## RECOMMENDED OPERATING PARAMETERS

<b>Nominal Retention Ratings</b>	0.2 to 0.8 micron
<b>Maximum Operating Pressure</b>	VR Cartridges : 2.4 bar maximum cartridge pressure drop. BioCap® 25 Capsules : 2.75 bar maximum inlet pressure, 2.4 bar maximum media pressure drop. BioCap® 1000 / 2000 Series Capsules : 5.5 bar at 25°C, 4.1 bar at 60°C, maximum inlet pressure ; 2.4 bar maximum cartridge pressure drop.
<b>Maximum Operating Temperature</b>	Zeta Plus® VR 02, 03, 05, 06 & 07 Cartridges : 82°C. Zeta Plus® VR 01 & 04 Cartridges : 60°C. BioCap® 25 Capsules : 40°C. BioCap® 1000 / 2000 Capsules : 60°C.
<b>Recommended Pre-use Rinse</b>	54 l/m <sup>2</sup> @ 20 lpm/m <sup>2</sup>
<b>Sterilization Parameters</b>	All media grades can be autoclaved or in-situ steam sterilized, except VR 01 & VR 04 (refer to Table 1.) Zeta Plus® VR Cartridges : In-situ steam sterilized or autoclaved for 30 minutes at 121° C, 1 cycle. BioCap® 25 : Autoclave for 30 minutes at 121° C (1 cycle). BioCap® 1000 / 2000 : Autoclave for 30 minutes at 121° C (up to 3 cycles).
<b>Operating FLUX range</b>	1.2 - 12 lpm/m <sup>2</sup>

## SURFACE AREA FOR ZETA PLUS® VR FILTER PRODUCTS :

Nominal Surface Area	
BioCap® 25 Capsule	25 cm <sup>2</sup>
BioCap® 1000 Capsule	650 cm <sup>2</sup>
BioCap® 2000 Capsule	1300 cm <sup>2</sup>
45109 (8" diameter cartridge, 8-cell)	2600 cm <sup>2</sup>
45167 (8" diameter cartridge, 7-cell, O-ring, plug-in)	2300 cm <sup>2</sup>
Z8FA2NPX2 (8" diameter cartridge, 2-cell plug-in)	650 cm <sup>2</sup>
Z8FA4NPX2 (8" diameter cartridge, 4-cell plug-in)	1300 cm <sup>2</sup>
45264 (12" diameter cartridge, 7-cell bodyfeed)	0,67 m <sup>2</sup>
45245 (12" diameter cartridge, 16-cell)	1.5 m <sup>2</sup>
Z16P (16", 14-cell, netted)	3.20 m <sup>2</sup>

## CARTRIDGE CONSTRUCTION :



### VR SERIES ORDERING GUIDE - DISPOSABLE CARTRIDGE AND CAPSULE PRODUCTS

#### 8" DIAMETER CARTRIDGES

CATALOGUE NUMBER	GASKET	VR SERIES
<b>45109</b> (8", 8-Cell)	11 - Nitrile	<b>VR 01 VR 05</b> <b>VR 02 VR 06</b> <b>VR 03 VR 07</b> <b>VR 04</b>
	13 - Fluorocarbon	
	14 - EPR	
	22 - Silicone	
<b>45167</b> (8", 7-Cell, O-ring Plug-in)	01 - Nitrile	<b>VR 01 VR 05</b> <b>VR 02 VR 06</b> <b>VR 03 VR 07</b> <b>VR 04</b>
	02 - EPR	
	03 - Fluorocarbon	
	04 - Silicone	
<b>Z8FA2NP</b> (8", 2-Cell, Plug-in)	A - Silicone	<b>VR 01 VR 05</b> <b>VR 02 VR 06</b> <b>VR 03 VR 07</b> <b>VR 04</b>
	B - Fluorocarbon	
<b>Z8FA4NP</b> (8", 4-Cell Plug-in)	C - EPR	
	D - Nitrile	
<b>PKG. CODE</b>		2

#### 16" DIAMETER CARTRIDGES

CATALOGUE NUMBER	GASKET	VR SERIES
<b>Z16P</b> (16", 14-Cell)	A - Silicone	<b>VR 01 VR 05</b> <b>VR 02 VR 06</b> <b>VR 03 VR 07</b> <b>VR 04</b>
	B - Fluorocarbon	
	C - EPR	
	D - Nitrile	

#### BIOCAP® CAPSULES

CATALOGUE NUMBER	VR SERIES
<b>BC 0025L</b> (Luer) - 25 cm <sup>2</sup> area	<b>VR 01 VR 05</b> <b>VR 02 VR 06</b> <b>VR 03 VR 07</b> <b>VR 04</b>
<b>BC 0025S</b> (Sanitary) - 25 cm <sup>2</sup> area	
<b>BC 1000A</b> (single pack) - 650 cm <sup>2</sup> area	
<b>BC 1000B</b> (3 pack) - 650 cm <sup>2</sup> area	
<b>BC 2000A</b> (single pack) - 1300 cm <sup>2</sup> area	
<b>BC 2000B</b> (3 pack) - 1300 cm <sup>2</sup> area	

#### 12" DIAMETER CARTRIDGES

CATALOGUE NUMBER	GEOMETRIC VARIATION	GASKET	VR SERIES
<b>45264</b> (12", 7-Cell Bodyfeed)	01 - Standard Polypropylene	A - Silicone	<b>VR 01 VR 05</b> <b>VR 02 VR 06</b> <b>VR 03 VR 07</b> <b>VR 04</b>
<b>45245</b> (12", 16-Cell)		B - Fluorocarbon	
	C - EPR		
	D - Nitrile		

#### Important Notice

CUNO MAKES NO WARRANTIES, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. Since a variety of factors can affect the use and performance of a CUNO product in a particular application, some of which are uniquely within the user's knowledge and control, user is responsible for determining whether or not the CUNO product is fit for a particular purpose and suitable for user's method of application.

#### Limitation of Remedies and Liability

If the CUNO product is proved to be defective, THE EXCLUSIVE REMEDY, AT CUNO'S OPTION, SHALL BE TO REFUND THE PURCHASE PRICE OR TO REPAIR OR REPLACE THE DEFECTIVE PRODUCT. CUNO shall not otherwise be liable for loss or damages, whether direct, indirect, special, incidental or consequential, regardless of the legal theory asserted, including, but not limited to, contract, negligence, warranty or strict liability.

#### WARRANTY

Seller warrants its equipment against defects in workmanship and material for a period of 12 months from date of shipment from the factory under normal use and service and otherwise when such equipment is used in accordance with instructions furnished by Seller and for purposes disclosed in writing at the time of purchase, if any. Any unauthorised alteration or modification of the equipment by Buyer will void this warranty. Seller's liability under this warranty shall be limited to the replacement or repair, F.O.B., point of manufacture, of any defective equipment or part which, having been returned to the factory, transportation charges prepaid, has been inspected and determined by Seller to be defective. THIS WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EITHER EXPRESSED OR IMPLIED, AS TO DESCRIPTION, QUALITY, MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE OR USE, OR ANY OTHER MATTER. Under no circumstances shall Seller be liable to Buyer or any third party for any loss of profits or other direct or indirect costs, expenses, losses or consequential damages arising out of or as a result of any defects in or failure of its products or any part or parts thereof or arising out of or as a result of parts or components incorporated in Seller's equipment but not supplied by the Seller.



#### 3M Europe CUNO Filtration

Hermeslaan 7  
1831 Diegem  
Belgium

Phone: +32-2-7224500

Fax: +32-2-7224518

E-mail: infocuno-europe@mmm.com

Web site address: www.3m.eu

For more contact addresses visit our website [www.cuno.com/international](http://www.cuno.com/international) or [www.3m.eu/filtration](http://www.3m.eu/filtration).

Data may be subject to change without further notice.

© 3M 2007. All rights reserved.

LITZPVR.EU